A1C

HbA1c analyzer

A Better A1C Test







Specifications

Testing Principle	Affinity Chromatography	Data Storage	10,000 test results
	, , , ,	NA/a: alab	.41/m
Testing Range	4%-16%	Weight	<1Kg
Precision	CV≤3%	Analyzer Dimension	200*160*80 mm
Blood Sample	Capillary/Venous Blood	Operating Condition	Temperature: 10°C-40°C; Humidity:30%-70%
Testing Time	< 3mins	Storage Condition	Temperature: -10°C-40°C; Humidity:<80%
Sample Volume	10µl	Test Kit	25 tests/pack; 2°C-30°C; 12months

Operation



Step 1: Adds 10µl R1-A to R1-B with capillary tube, mixes well



O 2mir

Step 2: Adds 10µl blood sample to the mixed solution, mixes well: incubates for 2mins



Step 3: Adds one drop of mixed solution and one drop of R2, waits 10s and read result

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Ultrasound Imaging
- In-Vitro Diagnostics

- Patient Monitoring
- Point-of-Care Testing
- Veterinary

• OB/GYN

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen 518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. \mid 9918 Via Pasar, San Diego, CA 92126

+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com