


PADECG

A Mobile ECG Solution

Versatile
Mobile ECG Solution




EDAN

A world of potential

PADECG A Mobile ECG Solution

PADECG consists of a professional ECG software and a wireless ECG acquisition module. It provides a mobile solution to meet the demands of caregivers in a mobile treatment environment. EDAN's reliable SEMIP algorithm is also provided in it, which allows you to have a fast and confident ECG assessment before sending the ECG report to Smart ECG Viewer for further analysis.

PADECG Paperless Workflow



- User-friendly Interfaces
- Reliable Interpretation
- Extensive Communication

Touch Screen

- Operating PADECG is just as simple as using your smart phone. Convenient gestures such as long pressing, un-pinching, and dragging are all supported to enable users a better experience.

ECG Sampling

- Both pre-sampling and real-time sampling modes are offered to cater to different situation, no matter it is an emergent case, general ward rounds, or home visiting.


Filter

- Armed with the comprehensive filters, you can acquire superb ECG traces with full confidence just within several touches.

Algorithm

- With EDAN's professional ECG algorithm SEMIP, PADECG gives accurate diagnosis results and offers doctors a reliable reference.

★ *Tablet is not included in PADECG. Customers need to purchase tablet by themselves.*



Global Headquarters:

Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen
518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:

EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126
+1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-ECG-PADECG-Android
V1.0-20200312